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### **English** Bipolar vessel sealing and coagulation instrument

#### **Attention**

Please read all information contained in this insert.  
Incorrect handling and care as well as misuse can lead to premature wear of surgical instruments.

#### **Intended Use**

The POWEREDGE bipolar vessel sealing and coagulation instrument with cutting mechanism is intended for use in minimally invasive and especially laparoscopic surgical procedures. The device is to be passed through a 5.5 mm laparoscopic cannula. The Instrument is designed to grasp, compress, seal or coagulate selected vessels and vascular bundles or tissue so that they can then be mechanically separated with the appropriate POWEREDGE blade.  
It must be connected to the bipolar output of an electrological generator. As indicated bipolar coagulation current may be selectively applied to the tissue.

**Maximum output voltage of the generator. U<sub>max</sub>:**  
400 V<sub>p</sub>

#### **Bissinger accessories:**

Bissinger POWEREDGE-blade REF 82710300  
Bissinger bipolar cable REF 80100xxx.

**Attention: Instruments for electrosurgery should be used only by persons who have been specially trained in the use of such instruments.**

#### **Contraindications**

**Incidents which have been reported in connection with the use of bipolar systems:**  
Unintended activation with resulting tissue injury on the wrong spot and/ or damage to the equipment.  
Fire in connection with surgical drapes and other inflammable materials.  
Alternating current paths leading to burns on spots where the patient or user comes into contact with components without insulation.  
Explosions caused by sparks in the proximity of inflammable gases.  
Perforation of organs. Sudden severe bleedings.

#### **Use and safety instructions**

- The non-observance of the present use and safety instructions may lead to injuries, malfunctions or other unexpected incidents.
- Before initial use and any other use, all instruments have to be completely cleaned, disinfected and sterilised and their functioning has to be checked.
- The POWEREDGE-blade REF 82710300 is intended for one-time use and must be disposed after disassembly.
- It is very important to check each surgical instrument for visible damage and wear, such as cracks, breaks or insulation defects before each use. In particular areas, such as blades, tips, notches, locking and blocking devices as well as all mobile parts and insulations have to be checked carefully.
- Never use damaged instruments.
- Never use the instruments in the presence of flammable or explosive substances.
- The instrument may not be laid down on the patient.
- Coagulation should only be performed if the contact surfaces are visible and ensure a good contact to the tissue selected for coagulation. Do not touch any other metallic instruments, trocar sleeves, optics or things like that during use.
- Only use the cutting mechanism when the enclosed tissue or blood vessel is completely coagulated or sealed in order to avoid bleeding and post-bleeding. Observe the use and safety instructions of the manufacturer of the high-frequency surgical device.

#### **Assembly and Operation**

Follow the pictogram PIC-827 for assembly and disassembly of the instrument.  
Once correctly assembled and connected to power, the device may be used in either the right or the left hand.  
To close jaws: compress (grip) handle.  
To open jaws: release (grip) handle.  
To activate the blade: depress the trigger.  
To deactivate the blade: release the trigger  
Coagulation current is activated by a footswitch that is part of your electrosurgical unit.

#### **Reprocessing**

Due to the product design, the raw materials used and the intended purpose it is not possible to determine a precise limit with regard to the maximum possible number of reprocessing cycles. The serviceable life of the instruments is determined by their function as well as by a careful handling.

The POWEREDGE-blade REF 82710300 is only intended for one-time use and must not be reused.  
Instruments for electrosurgery are by nature subject to increased wear depending on the type and time of use.

#### **Preparation and transport**

Disassemble the instrument following pictogram PIC 827 immediately after each use. Discard the POWEREDGE blade REF 82710300. The instrument must be cleaned and disinfected in a disassembled and open condition. Remove coarse dirt from the instruments. Do not use fixation agents or hot water (>40°C). Storage and transport of the instruments to the reprocessing location must be ensured in a sealed container.

#### **Manual pre-cleaning**

1. Immerse the instrument in cold water for 5 minutes.
2. Brush the instruments under cold water until all visible impurities are removed.
3. Flush each surface, which is not accessible with the brush, with a water jet pistol for 10 seconds (4 bar water pressure).
4. Place the instruments in an ultrasonic bath with a 0.5% alkaline-enzymatic cleaning detergent (MediClean forte, Dr. Weigert Hamburg). Ultrasonic must be applied for 15 minutes at 40°C/104°F. Make sure that the instruments are completely wet. The electrode must be in an opened condition when placing in the bath.
5. Remove the instrument and rinse completely with cold tap water to remove the cleaning detergent.

#### **Machine reprocessing**

##### **Cleaning**

Place the instruments in a basket on the insert module or on the inserts of the MIS module and start the cleaning process.

1. 4 min. pre-rinse with cold water.
2. Drain
3. 10 min. wash at 45°C with 0.5% of the alkaline-enzymatic cleaning detergent (Neodisher MediClean forte, Dr. Weigert Hamburg)
4. Drain
5. 3 min. neutralization with cold demineralized water.
6. Drain
7. 2 min. post-rinse with cold demineralized water.
8. Drain

##### **Disinfection**

Machine operated thermal disinfection has to be carried out in consideration of the national requirements with regard to the A0 value (see ISO 15883).

##### **Drying**

Dry the outside of the instruments by carrying out a drying cycle of the cleaning/disinfection machine.  
If necessary, manual drying may additionally be carried out using a lintfree cloth. Dry cavities by blowing with sterile compressed air.

##### **Manual reprocessing**

Cannot be applied to this instrument.

##### **Functional test and packaging**

Perform a visual inspection for cleanliness and good condition of the material. If necessary, repeat the reprocessing process until the instrument is optically clean.  
Use harmless medical white oil (e.g. paraffin oil acc. to DAB or Ph.Eur. or USP) to care for the thread of the instruments after cleaning and before assembling.  
Assemble the instrument following PIC 827, which will require the use of an additional POWEREDGE blade REF 82710300.  
Check the functionality of the instrument based on the application notes. Packaging has to comply with ISO 11607 and EN 868 standards for the packaging for terminally sterilized medical devices.

##### **Sterilisation**

Sterilisation of the products with fractional pre-vacuum procedure (in accordance with ISO 13060 / ISO 17665) in consideration of the respective national requirements.  
- 3 pre-vacuum phases with a pressure of at least 60 mbar.  
- Heating up to a sterilisation temperature of min. 132°C and max. 137°C.  
- Shortest exposure time: 3 min.  
- Drying time: at least 10 min.

##### **Storage**

Sterilised instruments have to be stored in a dry, clean and dust-free area at moderate temperatures from 5°C to 40°C.

## Repairs

Never attempt to perform repairs yourself. Service and repair work may only be performed by persons qualified and trained accordingly. For any question on these matters, please contact either the manufacturer or your medico-technical department.

**Attention:** Defect products must pass the complete reprocessing process before being returned for repair.

## Information on the validation of the reconditioning

The following testing instructions, materials and equipment have been used for validation:

### Cleaning agents (for machine use and ultrasonic):

- Neodisher MediClean forte; Dr. Weigert Hamburg
- Cleaning and disinfection device:
- Miele G7836 CD
- MIS rack E450
- Steam steriliser:
- Selectomat HP 666-1HR

For details, see report.

SMP GmbH # 20214-2 (machine cleaning)  
SMP GmbH # 25714-2 (sterilisation)

If the chemicals and machines described above are not available, the user has to validate the used process accordingly.

## Handling

During transport, cleaning, care, sterilisation and storage, all surgical instruments should be handled with maximum care.

This applies particularly to blades, fine tips and other sensitive areas.

## Warranty

Günter Bissinger Medizintechnik GmbH exclusively supplies tested and faultless products to the customers. All products are designed and manufactured to comply with maximum quality requirements. We refuse any liability for products which, compared to the original product, have been modified, misused or handled or used in an inexpert way.

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SMP GmbH # 20214-2 (limpieza mecánica)
SMP GmbH # 25714-2 (esterilización)
DQS Medizinprodukte GmbH August-Schanz-Straße 21 60433 Frankfurt, Germany