

bissinger

PLASMA LOOP

TRUE BIPOLAR RESECTOSCOPE

			
<p>Art. No. 86700002 5</p>	<p>Art. No. 86700004 5</p>	<p>Art. No. 86700014 5</p>	<p>Art. No. 86700006 5</p>
<p>PLASMALOOP Bip. Elektrode Schlinge 90° gewinkelt PLASMALOOP Bip. electrode Loop 90° angled</p>	<p>PLASMALOOP Bip. Elektrode Schlinge 60° gewinkelt PLASMALOOP Bip. electrode Loop 60° angled</p>	<p>PLASMALOOP Bip. Elektrode Kugel PLASMALOOP Bip. electrode ball</p>	<p>PLASMALOOP Bip. Elektrode Schlinge 0° gerade PLASMALOOP Bip. electrode Loop 0° straight</p>
<p>NEW</p> 	<p>NEW</p> 		
<p>Art. No. 86700008 5</p>	<p>Art. No. 86700010 5</p>	<p>Art. No. 86700012 5</p>	<p>Art. No. 86700016 5</p>
<p>PLASMALOOP Bip. Elektrode Big Loop 90° PLASMALOOP Bip. electrode Big Loop 90°</p>	<p>PLASMALOOP Bip. Elektrode Big Loop 60° PLASMALOOP Bip. electrode Big Loop 60°</p>	<p>PLASMALOOP Bip. Elektrode Messer PLASMALOOP Bip. electrode knife</p>	<p>PLASMALOOP Bip. Elektrode Walze PLASMALOOP Bip. electrode roller</p>

PLASMALOOP

TRUE BIPOLAR RESECTOSCOPE

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In general, transurethral resection (TUR) requires high power in order to ensure a constant cutting quality.

In monopolar procedures the high current flow is passed from the surgical field through surrounding tissue and the body of the patient back to the HF surgical unit via a neutral electrode, whereas in bipolar procedures current flow is limited to the operating field only.

The bipolar Plasmaloop system greatly reduces the risks associated with the use of conventional monopolar systems, risks such as uncontrolled current flow are eliminated. With Plasmaloop, current flow is limited to the space between the two distally placed electrodes and occurs only in the area visible to the surgeon.

As bipolar resection is performed in saline solution, the risk of a TUR syndrome occurrence is practically eliminated.

The Plasmaloop bipolar electrode arrangement ensures outstanding ignition performance, with precise, reliable and consistent cutting and coagulation quality throughout the whole procedure.

This specially designed arrangement of the two distal electrodes also ensures that there is no current flow along the sheath, thus preventing unwanted thermal tissue damage caused by leakage current or damaged insulation.

The Plasmaloop system design allows simple insertion and removal of the electrode.

The ergonomic handle, the simple and self-explanatory detachability of the instrument and the cable connector being positioned parallel to the working axis enable safe and efficient working.

All components of the Plasmaloop system are made from high quality materials that easily withstand current cleaning and sterilisation requirements.

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Medizintechnik GmbH**