




Xclusive

Bipolare Pinzette Claris Non-Stick
Bipolar forceps Claris Non-Stick



Xclusive

CLARIS Non-Stick



-Neuartiges Design für präzise Anwendungen-

Die extrem schlanken Pinzettenschenkel mit Antihaft-Spitzen aus echtem Silber, erlauben das Arbeiten in schwer zugänglichen Bereichen ebenso, wie eine präzise Anwendung im fragilen Gewebe.

Das optimal gewählte Verhältnis von Schenkel zu Spitze erlaubt dem Operateur uneingeschränkte Sicht auf die Silberspitzen. Dadurch wird im Zusammenspiel mit der integrierten Schenkelführung ein zielgenaues koagulieren und präparieren gewährleistet.

-New design for precise applications-

The ultra-slim forceps shanks with non-stick tips, made of silver, enhances working in hard to reach areas as well as ensuring precise application in fragile tissue.

The optimally selected arm design and angle provides the surgeon with an unrestricted view of the silver tips. This ensures precise preparation and coagulation when combined with the integrated leg guide.



“ We believe that Bissinger CLARIS Xclusive non-stick-forceps are the most ideal instruments for precise operations. ”

- Dr. Haishi Zhang -
(Huashan Hospital, Fudan University, Shanghai)

Spitzenbreite
Width of tips

Anschluß
Plug

0,3 mm x 3 mm



240 mm

80524337



**Günter Bissinger
Medizintechnik GmbH**

Hans-Theisen-Str.1
79331 Teningen/Germany

Tel. +49 7641 9 14 33 0

Fax +49 7641 9 14 33 33

eMail: info@bissinger.com

Internet: <http://www.bissinger.com>