


80522348

CLARIS Xclusive Non-Stick Bipolar Forceps, bayonet



## Technical specifications

### Dimensions and Weights

|                        |        |
|------------------------|--------|
| Total length           | 220 mm |
| Total length in inches | 8,66 " |
| Tip width              | 0,3 mm |

### Basic Specifications

|               |           |
|---------------|-----------|
| Forceps shape | bayonet   |
| Tip shape     | straight  |
| Connector     | Flat Plug |

### Cleaning and Reprocessing

|                         |                     |
|-------------------------|---------------------|
| Reusable                | Yes                 |
| Manual reprocessing     | Yes, related to IFU |
| Washing machine         | 95 °C               |
| Autoclave sterilization | max. 137 °C         |

### Approvals

|                          |           |
|--------------------------|-----------|
| Directive 93/42 EEC      | Class IIb |
| Regulation (EU) 2017/745 | N/A       |
| FDA Approval             | Class II  |

### Additional Product Information

|                           |                |
|---------------------------|----------------|
| Accessories rated voltage | max. 250 Vp    |
| Packaging unit            | 1              |
| GTIN 14                   | 04250418749039 |