


30-0564-00

Cystoscope



## Technical specifications

### Dimensions and Weights

### Cleaning and Reprocessing

|                         |                     |
|-------------------------|---------------------|
| Reusable                | Yes                 |
| Manual reprocessing     | Yes, related to IFU |
| Washing machine         | 95 °C               |
| Autoclave sterilization | max. 137 °C         |

### Approvals

|                          |         |
|--------------------------|---------|
| Directive 93/42 EEC      | Class I |
| Regulation (EU) 2017/745 | N/A     |
| FDA Approval             | N/A     |

### Additional Product Information

|                |                |
|----------------|----------------|
| Packaging unit | 1              |
| GTIN 14        | 04250418740821 |