

82460326

## MITHRAS Bipolar Coagulation Forceps, Ring Handle



### Technical specifications

#### Dimensions and Weights

Total length	150 mm
Diameter	3 mm

#### Basic Specifications

Electrode shape	Scissor, curved
-----------------	-----------------

#### Cleaning and Reprocessing

Reusable	Yes
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

#### Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	Class II

#### Additional Product Information

Accessories rated voltage	max. 300 Vp
Packaging unit	1
GTIN 14	04250418748346