

89217024

Monopolar Electrode Handle



Technical specifications

Dimensions and Weights

| | |
|------------|----|
| Dimensions | mm |
|------------|----|

Basic Specifications

| | |
|----------------|--------|
| Shaft diameter | 2,4 mm |
|----------------|--------|

Cleaning and Reprocessing

| | |
|----------|-----|
| Reusable | Yes |
|----------|-----|

| | |
|---------------------|---------------------|
| Manual reprocessing | Yes, related to IFU |
|---------------------|---------------------|

| | |
|-----------------|-------|
| Washing machine | 95 °C |
|-----------------|-------|

| | |
|-------------------------|-------------|
| Autoclave sterilization | max. 137 °C |
|-------------------------|-------------|

Approvals

| | |
|---------------------|-----------|
| Directive 93/42 EEC | Class IIb |
|---------------------|-----------|

| | |
|--------------------------|-----|
| Regulation (EU) 2017/745 | N/A |
|--------------------------|-----|

| | |
|--------------|-----|
| FDA Approval | N/A |
|--------------|-----|

Additional Product Information

| | |
|---------------------------|------------|
| Accessories rated voltage | max. 5 kVp |
|---------------------------|------------|

| | |
|----------------|---|
| Packaging unit | 1 |
|----------------|---|

| | |
|---------|----------------|
| GTIN 14 | 04250418728164 |
|---------|----------------|