

85904001

Monopolar Hook for MIS



Technical specifications

Dimensions and Weights

| | |
|--------------------|--------|
| Total length | 290 mm |
| Total length in cm | 29 cm |
| Diameter | 3 mm |
| Angle | 0 ° |

Basic Specifications

| | |
|-----------------|--------|
| Electrode shape | L-Hook |
|-----------------|--------|

Cleaning and Reprocessing

| | |
|-------------------------|---------------------|
| Reusable | Yes |
| Manual reprocessing | Yes, related to IFU |
| Washing machine | 95 °C |
| Autoclave sterilization | max. 137 °C |

Approvals

| | |
|--------------------------|-----------|
| Directive 93/42 EEC | Class IIb |
| Regulation (EU) 2017/745 | N/A |
| FDA Approval | N/A |

Additional Product Information

| | |
|---------------------------|----------------|
| Accessories rated voltage | max. 2 kVp |
| Packaging unit | 1 |
| GTIN 14 | 04250418742566 |