


**84650413**

Needle Holder Ø 5 mm, Curved Handle



## Technical specifications

### Dimensions and Weights

Total length	330 mm
Diameter	5 mm

### Basic Specifications

#### Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

#### Approvals

Directive 93/42 EEC	Class I
Regulation (EU) 2017/745	N/A
FDA Approval	N/A

#### Additional Product Information

Packaging unit	1
GTIN 14	04250418749145