

**82400111**

**POWERGRIP / ORBITARIS Bipolar Coagulation Forceps, Electrode**



## Technical specifications

### Dimensions and Weights

Total length	340 mm
Diameter	5 mm

### Basic Specifications

Electrode shape	Grasping Forceps with 2x3 teeth
-----------------	---------------------------------

### Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

### Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	Class II

### Additional Product Information

Accessories rated voltage	max. 500 Vp
Packaging unit	1
GTIN 14	04250418741446