


**85530004**

**SLIMLINE Bipolar Coagulation Forceps, Shaft**



## Technical specifications

### Dimensions and Weights

|              |        |
|--------------|--------|
| Total length | 340 mm |
| Diameter     | 3,4 mm |

### Cleaning and Reprocessing

|                         |                     |
|-------------------------|---------------------|
| Reusable                | Yes                 |
| Manual reprocessing     | Yes, related to IFU |
| Washing machine         | 95 °C               |
| Autoclave sterilization | max. 137 °C         |

### Approvals

|                          |           |
|--------------------------|-----------|
| Directive 93/42 EEC      | Class IIb |
| Regulation (EU) 2017/745 | N/A       |
| FDA Approval             | Class II  |

### Additional Product Information

|                |                |
|----------------|----------------|
| Packaging unit | 1              |
| GTIN 14        | 04250418723688 |